

EDAN

C3/C6

Video Colposcope



EDAN Edan Instruments, Inc.

• 3/F – B, Nanshan Medical Equipments Park, Nanhai Rd 1019#, Shekou, Nanshan Shenzhen, 518067 P.R. China
• Tel: +86-755-26898326 • Fax: +86-755-26898330 • www.edan.com.cn • Email: info@edan.com.cn

• All rights reserved. Features and specifications are subject to change without notice



ENG-OB & GYN-C3&C6-V1.2-20140415

Care for Health

With its bright LED cold lighting system, high-resolution camera and electronic green filter, C3/C6 Video Colposcope realizes full-screen display of precise cervical examination, offering crucial and accurate information for clinical decision making to the gynecologist.

Features

- LED cold lighting system reproduces an original view of tissue
- Advanced and fast auto-focusing system ensures clear image
- Remote capture control and user-friendly keyboard layout for easy operation
- Electronic green filter enhances vascular image display without compromise of illumination
- Acetic acid and iodine reaction test timer to capture the right and precise image
- Real-time display of magnification
- Powerful image display and data management software

R-way™

The R-way function observes the variation of "color" on cervix to interpret the performance of cervical diseases

- Saline application – 'R for Red'**
To observe the area with rich blood circulation, this reflects in red.

Acetic acid 5% – 'W for White'
To observe increase in nucleus amount, this reflects in white.

Green filter – 'A for Atypical'
To better observe atypical vascular, this defines as A.

Lugol's iodine solution – 'Y for Yellow'
To observe decrease in cell plasma, this reflects in yellow.
- CIN2+**



C3/C6 Video Colposcope

We also offer you patient management software especially designed for the needs of gynecologists, which enables you to take full advantage of the Video Colposcope. The software enables you to capture, review, analyze, store, print, and export colposcopic images quickly and easily.

- Patient tracking and data management
- Built-in marking system for analysis
- Export report in PDF format
- Print report on A4 size paper via external printer

